



Eldex® Eldex Laboratories, Inc.
 30 Executive Court
 Napa, CA 94558-6278
 USA

Tel: (707) 224-8800
 (800) 969-3533
 Fax: (707) 224-0688
 email: service@eldex.com

Instructions

MicroPro Quaternary/Ternary Connections

Overview

A ternary or quaternary MicroPro pumping system is based on expansion of a binary gradient system. The binary gradient system is a dual syringe pump system in a single chassis. Additional solvents are added to the binary system by adding additional single syringe “slave” units.

Figure 1 illustrates the intended configuration of a ternary system. A quaternary system differs only in the addition of another single syringe slave unit, to the left of the ternary pump.




Figure 1: Ternary Gradient Illustration

Plumbing

Systems are provided with the tubing required to connect from port 6 on the high pressure valve of the slave pump to the mixer port on the binary gradient system. (see Figure 2).




Figure 2: Plumbing Connections

Electrical

Address Port Settings

Figure 3 shows the location of the address port. A binary gradient system is set to position 2. The ternary slave pump is set to position 5. The quaternary slave pump is set to position 6.




Figure 3: Address Port/RS485

the binary gradient system and the slave units.

Additional Notes

Consult the MicroPro User Manual for additional details about voltage settings, other electrical and plumbing connections. Note that the I/O and Analog I/O boards on slave units perform no function.

Power

The power to the slave pumps should be turned on before the binary gradient system power is turned on.

Communication Cable

A standard, nine-pin cable is provided with a ternary or quaternary slave pump. Connection is made between the RS485 ports (see Figure 3) on